

## Product datasheet

# Anti-Fascin-3 antibody ab91484

## 2 图像

### 概述

产品名称	Anti-Fascin-3抗体
描述	兔多克隆抗体to Fascin-3
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide selected from the N-terminal region of Human Fascin-3 conjugated to KLH (NP_065102).
阳性对照	NCI-H460 cell line lysate and mouse testis tissue lysate.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	ab91484 is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab91484** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
ELISA		
应用说明	ELISA: 1/1000.	

WB: 1/50 - 1/100. Predicted molecular weight: 57 kDa.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标


功能	Acts as an actin bundling protein.
组织特异性	Expressed in testis.
序列相似性	Belongs to the fascin family.
细胞定位	Cytoplasm > cytoskeleton.

## 图片



Anti-Fascin-3 antibody (ab91484) at 1/50 dilution + NCI-H460 cell line lysate at 35  $\mu$ g

**Predicted band size : 57 kDa**



Anti-Fascin-3 antibody (ab91484) at 1/50 dilution + Mouse testis tissue lysate at 35  $\mu$ g

**Predicted band size : 57 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors