

Product datasheet

Anti-Fascin-3 antibody ab91484

2 图像

概述

产品名称	Anti-Fascin-3抗体
描述	兔多克隆抗体to Fascin-3
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide selected from the N-terminal region of Human Fascin-3 conjugated to KLH (NP_065102).
阳性对照	NCI-H460 cell line lysate and mouse testis tissue lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	ab91484 is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab91484** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
ELISA		

应用说明

ELISA: 1/1000.

WB: 1/50 - 1/100. Predicted molecular weight: 57 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


功能 Acts as an actin bundling protein.

组织特异性 Expressed in testis.

序列相似性 Belongs to the fascin family.


细胞定位 Cytoplasm > cytoskeleton.

图片



Anti-Fascin-3 antibody (ab91484) at 1/50 dilution + NCI-H460 cell line lysate at 35 μ g

Predicted band size: 57 kDa



Anti-Fascin-3 antibody (ab91484) at 1/50 dilution + Mouse testis tissue lysate at 35 μ g

Predicted band size: 57 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors