

Product datasheet

Anti-GC1q R antibody - Azide free ab37851

1 图像

概述

产品名称	Anti-GC1q R抗体- Azide free
描述	鸡多克隆抗体to GC1q R - Azide free
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat, Dog, Human
免疫原	Synthetic peptide, corresponding to C terminal amino acids 74-282 of Human GC1q R

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: None Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgY

应用

Our [Abpromise guarantee](#) covers the use of **ab37851** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 31 kDa.

靶标

功能	Binds to the globular "heads" of C1Q thus inhibiting C1 activation.
序列相似性	Belongs to the MAM33 family.
细胞定位	Mitochondrion matrix. Nucleus. Might also be nuclear in some cell types.

## Anti-GC1q R antibody - Azide free 图像



Anti-GC1q R antibody - Azide free (ab37851)  
at 1/2000 dilution + E. coli-derived fusion  
protein at 2.5 µg

**Predicted band size :** 31 kDa

**Observed band size :** 80 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors