

Product datasheet

Anti-GM-CSF antibody [59] ab649

1 References

概述

| | |
|--------------|--|
| 产品名称 | Anti-GM-CSF抗体[59] |
| 描述 | 小鼠单克隆抗体[59] to GM-CSF |
| 宿主 | Mouse |
| 特异性 | This antibody reacts with natural and recombinant hGM-CSF. It shows no cross-reaction with other human cytokines or growth factors tested such as IL-1 beta, IL-8, IL-16, TGF-beta 1, GM-CSFR, TNF-alpha, and M-CSF. |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Purified recombinant hGM-CSF (human granulocyte macrophage colony stimulating factor) |
| 阳性对照 | <div style="border: 1px solid #ccc; padding: 5px; display: inline-block;"> <p>购买相配的WB阳性对照 Recombinant human GM-CSF protein ></p> </div> |
| 常规说明 | Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. If you would like help finding an alternative product, please do not hesitate to contact our scientific support team. |

性能

| | |
|-------------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 纯化说明 | Affinity purified from ascitic fluid. |
| 克隆 | 单克隆 |
| 克隆编号 | 59 |
| 骨髓瘤 | Sp2/0 |
| 同种型 | IgG1 |
| 轻链类型 | kappa |

应用

Our [Abpromise guarantee](#) covers the use of **ab649** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | Use a concentration of 1 µg/ml. At this concentration this antibody will visualize 0.2 µg/lane of hGM-CSF |

靶标

| | |
|-------|--|
| 功能 | Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes. |
| 序列相似性 | Belongs to the GM-CSF family. |
| 细胞定位 | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors