

## Product datasheet

# Anti-GPD2 antibody ab113058

★☆☆☆☆ 1 Abreviews 1 图像

### 概述

产品名称	Anti-GPD2抗体
描述	兔多克隆抗体to GPD2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 605-634 of Human GPD2.
阳性对照	ZR-75-1 cell line lysates

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab113058 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab113058** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB	★☆☆☆☆	1/100 - 1/500. Predicted molecular weight: 67 kDa.

### 图标


---

<b>通路</b>	Polyol metabolism; glycerol degradation via glycerol kinase pathway; glycerone phosphate from sn-glycerol 3-phosphate (anaerobic route): step 1/1.
<b>序列相似性</b>	Belongs to the FAD-dependent glycerol-3-phosphate dehydrogenase family. Contains 2 EF-hand domains.
<b>细胞定位</b>	Mitochondrion.

---

## 图片

---



Anti-GPD2 antibody (ab113058) at 1/100 dilution + ZR-75-1 cell line lysates at 35 µg

**Predicted band size:** 67 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors