

Product datasheet

Anti-GST antibody (FITC) ab3445

概述

| | |
|-------|--|
| 产品名称 | Anti-GST抗体(FITC) |
| 描述 | 兔多克隆抗体to GST (FITC) |
| 偶联物 | FITC. Ex: 493nm, Em: 528nm |
| 经测试应用 | 适用于: ELISA, ICC/IF, WB |
| 免疫原 | Full length protein corresponding to Schistosoma japonicum GST. Database link: P08515 |
| 常规说明 | Molar F/P ratio is 6.1. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | PBS with 0.2% BSA, 0.1% sodium azide |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | Immunoaffinity purified using GST immobilized on a solid support. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab3445** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|--------|------|--|
| ELISA | | Use at an assay dependent dilution. |
| ICC/IF | | Use at an assay dependent dilution. |
| WB | | Use at an assay dependent dilution. Predicted molecular weight: 1 kDa. |

靶标

相关性

GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors