

## Product datasheet

# Anti-GST antibody (FITC) ab3445

### 概述

产品名称	Anti-GST抗体(FITC)
描述	兔多克隆抗体to GST (FITC)
宿主	Rabbit
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: ELISA, ICC/IF, WB
免疫原	Full length protein corresponding to Schistosoma japonicum GST. Database link: <a href="#">P08515</a>
常规说明	Molar F/P ratio is 6.1.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	PBS with 0.2% BSA, 0.1% sodium azide
纯度	Immunogen affinity purified
纯化说明	Immunoaffinity purified using GST immobilized on a solid support.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab3445** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 1 kDa.

## 靶标

---

**相关性** GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors