

Anti-Gemin 4 antibody ab92535

3 图像

概述

产品名称	Anti-Gemin 4抗体
描述	兔多克隆抗体to Gemin 4
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Chimpanzee, Gorilla, Orangutan 
免疫原	Synthetic peptide, corresponding to a region between residue 1010 and the C-terminus (residue 1058) fo Human Gemin 4 (NP_056536.1).
阳性对照	Human breast carcinoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: Tris buffered saline, 0.1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab92535于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		

应用说明 IHC-P: 1/100 - 1/500.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


相关性

GEMIN 4 is part of the SMN complex and plays an essential role in spliceosomal snRNP assembly in the cytoplasm. It is required for pre mRNA splicing in the nucleus.

细胞定位


Cytoplasmic and nuclear.

图片




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Gemin 4 antibody (ab92535)

ab92535 at 1/250 dilution staining Gemin 4 in Human breast carcinoma by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue. Detection: DAB staining.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Gemin 4 antibody (ab92535)

ab92535 at 1/250 dilution staining Gemin 4 in human colon carcinoma by Immunohistochemistry Formalin fixed paraffin-embedded tissue. Detection DAB staining.



ab92535 at 1/250 dilution staining Gemin 4 in Mouse renal cell carcinoma by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue. Detection: DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Gemin 4 antibody (ab92535)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors