


Product datasheet

Anti-Gemin 4 antibody ab92535

3 图像

概述

产品名称	Anti-Gemin 4抗体
描述	兔多克隆抗体to Gemin 4
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Chimpanzee, Gorilla, Orangutan 
免疫原	Synthetic peptide, corresponding to a region between residue 1010 and the C-terminus (residue 1058) fo Human Gemin 4 (NP_056536.1).
阳性对照	Human breast carcinoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.1% BSA, Tris buffered saline
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab92535** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

IHC-P

应用说明 IHC-P: 1/100 - 1/500.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


相关性

GEMIN 4 is part of the SMN complex and plays an essential role in spliceosomal snRNP assembly in the cytoplasm. It is required for pre mRNA splicing in the nucleus.

细胞定位


Cytoplasmic and nuclear.

图片




ab92535 at 1/250 dilution staining Gemin 4 in Human breast carcinoma by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue. Detection: DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Gemin 4 antibody (ab92535)



ab92535 at 1/250 dilution staining Gemin 4 in human colon carcinoma by Immunohistochemistry Formalin fixed paraffin-embedded tissue. Detection DAB staining.

- Anti-Gemin 4 antibody (ab92535)



ab92535 at 1/250 dilution staining Gemin 4 in Mouse renal cell carcinoma by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue. Detection: DAB staining.

- Gemin 4 antibody (ab92535)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors