

Product datasheet

Goat Anti-Armenian Hamster IgG H&L (Biotin) ab5744

5 References 1 图像

概述

产品名称	山羊Anti-Armenian仓鼠IgG H&L (Biotin)
靶标种属	Armenian hamster
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Goat Serum, Armenian Hamster IgG and Armenian Hamster Serum. Cross reactivity against mouse has not been determined.
经测试应用	<b>适用于:</b> WB, ICC/IF, Dot blot, ELISA, Immunomicroscopy, IHC-P
免疫原	Armenian Hamster IgG whole molecule
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Goat IgG molecule.
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab5744** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000 - 1/10000.

ICC/IF	1/1000 - 1/5000.
Dot blot	Use at an assay dependent dilution.
ELISA	1/250000.
Immunomicroscopy	Use at an assay dependent dilution.
IHC-P	Use at an assay dependent dilution.

### Goat Anti-Armenian Hamster IgG H&L (Biotin) 图像



Ab5744 was used at dilution 1/500 with the primary antibody [ab33483](#) in IHC-P. See the review on [ab33483](#).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat polyclonal Secondary Antibody to Armenian Hamster IgG - H&L (Biotin) (ab5744)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors