

Product datasheet

Goat Anti-Chicken IgY H&L (10nm Gold) ab41511

1 References

概述

产品名称	山羊Anti-Chicken IgY H&L (10nm Gold)
靶标种属	Chicken
经测试应用	适用于: Electron Microscopy
免疫原	Chicken IgG whole molecule.
偶联物	Gold 10nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 15mM Sodium Azide Constituents: 1% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	This product has been immuno-affinity purified and immuno cross-absorbed to reduce non-specific reactions. The activity of each lot is determined using a dot-spot test system [Moeremans et al., J. Immunol., Methods 74, (1984), 353.]
克隆	多克隆
同种型	IgG
常规说明	This product is an immunogold reagent built around colloidal gold particles of 10nm. The particle population is monodisperse and thus shows minimal size variation and overlap. Typically, the coefficient of variance for the 10nm particle size conjugate shows less than 10% variation.

应用

Our [Abpromise guarantee](#) covers the use of **ab41511** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Electron Microscopy

应用说明

Electron Microscopy: Use at an assay dependent dilution. The method of production of this

product allows the small particle size of 10nm to be used instead of 5nm with the same sensitivity but increased visibility in the standard transmission electron microscope (See also ven Bergen en Henegouwen and Leunissen, Histochemistry, 85,(1986), 81.)

Light Microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors