

## Product datasheet

# Goat Anti-Donkey IgG H&L (FITC) ab6984

### 概述

产品名称	山羊抗驴IgG H&L (FITC)
宿主	Goat
靶标种属	Donkey
经测试应用	适用于: ICC/IF, IHC-Fr, Immunomicroscopy, Flow Cyt, ELISA, IHC-P
免疫原	Full length native Donkey IgG (purified).
偶联物	FITC. Ex: 493nm, Em: 528nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: 10mg/ml BSA, 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
纯度	Affinity purified
纯化说明	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons). Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Rabbit IgG.
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using donkey IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
克隆	多克隆

### 应用

Our [Abpromise guarantee](#) covers the use of **ab6984** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		
IHC-Fr		

应用	Ab评论	说明
Immunomicroscopy		
Flow Cyt		
ELISA		
IHC-P		

#### 应用说明

Flow Cyt: Use at an assay dependant dilution.  
 IHC-P: Use at an assay dependant dilution.  
 IHC-Fr: Use at an assay dependant dilution.  
 IM: Use at an assay dependant dilution.

Not tested in other applications.  
 Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors