

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG - (Fab')₂ (DyLight® 550), pre-adsorbed ab98758

1 图像

概述

产品名称	山羊F(ab') ₂ 抗小鼠IgG - (Fab') ₂ (DyLight® 550), 预吸附二抗	
宿主	Goat	
靶标种属	Mouse	
经测试应用	适用于: Flow Cyt, ICC/IF, IHC-P	
最小交叉反应	Human, Rat	more details
偶联物	DyLight® 550. Ex: 562nm, Em: 576nm	

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 550.
克隆	多克隆
同种型	IgG
常规说明	DyLight® 550 replaces DyLight® 549.

应用


Our [Abpromise guarantee](#) covers the use of **ab98758** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		1/50 - 1/200.

应用	Ab评论	说明
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

DyLight®-Goat F(ab')₂ polyclonal Secondary Antibody to Mouse IgG - (Fab')₂ (DyLight® 550), pre-adsorbed(ab98758)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors