

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG Fc (SureLight® Allophycocyanin) ab99903

概述

产品名称	山羊F(ab') ₂ 抗小鼠IgG Fc (SureLight® Allophycocyanin荧光)
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
偶联物	SureLight® Allophycocyanin. Ex: 652nm, Em: 657nm

性能

形式	Lyophilised: Lyophilized material is stable for up to one year. Reconstitute with 0.5 ml of ddH ₂ O. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 10mM Tris, 1 ug/ml pentachlorophenol pH 8.2
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab99903** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors