

## Product datasheet

# Goat F(ab')<sub>2</sub> Anti-Mouse IgG+IgM+IgA H&L (FITC) preadsorbed ab5999

### 3 References

#### 概述

<b>产品名称</b>	山羊F(ab') <sub>2</sub> 抗小鼠IgG+IgM+IgA H&L (FITC)预吸附二抗
<b>宿主</b>	Goat
<b>靶标种属</b>	Mouse
<b>特异性</b>	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein and anti-Goat Serum. No reaction was observed against anti-Goat IgG F(c), anti-Pepsin or Human Serum Proteins.
<b>经测试应用</b>	<b>适用于:</b> Flow Cyt, ICC/IF, ELISA
<b>最小交叉反应</b>	Human <a href="#">more details</a>
<b>免疫原</b>	Mouse IgG, IgA, and IgM whole molecules
<b>偶联物</b>	FITC. Ex: 493nm, Em: 528nm

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.01% Sodium Azide Constituents: 10mg/ml BSA (IgG and protease free), 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
<b>纯化说明</b>	This product was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
<b>共轭说明</b>	Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 4.0 moles FITC per mole of Goat IgG F(ab') <sub>2</sub>
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab5999** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		1/500 - 1/2500.
ICC/IF		1/1000 - 1/5000.
ELISA		1/10000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors