

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG+IgM+IgA H&L preadsorbed ab5998

概述

产品名称	山羊F(ab') ₂ 抗小鼠IgG+IgM+IgA H&L预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum. No reaction was observed against anti-Goat IgG F(c), anti-Pepsin or Human Serum Proteins. This product is suitable for the detection of all Mouse immunoglobulin classes, isotypes and chain combinations.
最小交叉反应	Human more details
免疫原	Mouse IgG, IgA, and IgM whole molecules

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 8.00 Preservative: 0.01% Sodium azide Constituents: 4.7625% Sodium borate, 0.146% EDTA, 0.435% Sodium chloride
纯化说明	This product was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
克隆	多克隆
同种型	IgG

应用

应用说明	Suitable for highly specific immunological methods requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
-------------	--

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors