

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG+IgM+IgA H&L preadsorbed ab5998

概述

| | |
|---------------|--|
| 产品名称 | 山羊F(ab') ₂ Anti-Mouse IgG+IgM+IgA H&L preadsorbed |
| 靶标种属 | Mouse |
| 特异性 | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum. No reaction was observed against anti-Goat IgG F(c), anti-Pepsin or Human Serum Proteins. This product is suitable for the detection of all Mouse immunoglobulin classes, isotypes and chain combinations. |
| 最小交叉反应 | Human more details |
| 免疫原 | Mouse IgG, IgA, and IgM whole molecules |

性能

| | |
|-------------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | Preservative: 0.01% Sodium Azide Constituents: 0.125M Sodium Borate, 0.075M Sodium Chloride, 0.005M EDTA. pH 8 |
| 纯化说明 | This product was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

| | |
|-------------|--|
| 应用说明 | Suitable for highly specific immunological methods requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user. |
|-------------|--|

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors