

Product datasheet

Goat F(ab')₂ Anti-Rabbit IgG - H&L, pre-adsorbed ab98504

概述

| | |
|--------|--|
| 产品名称 | 山羊F(ab') ₂ 抗兔IgG - H&L,预吸附二抗 |
| 宿主 | Goat |
| 靶标种属 | Rabbit |
| 经测试应用 | 适用于: ICC, IHC-P, ELISA, WB |
| 最小交叉反应 | Chicken, Cow, Goat, Horse, Human, Mouse, Rat, Sheep more details |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | Preservative: 0.09% Sodium Azide Constituents: PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab98504** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|---|
| ICC | | Use at an assay dependent dilution. |
| IHC-P | | 1/200 - 1/2000. |
| ELISA | | 1/100 - 1/500. Coating. ELISA: 1/1000 - 1/30000. Primary. |

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | 1/1000 - 1/10000. Colorimetric. WB: 1/1000 - 1/30000. Chemiluminescent. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors