

Product datasheet

Goat Anti-Human IgA H&L (FITC) ab8500

概述

产品名称	山羊抗人IgA H&L (FITC)
宿主	Goat
靶标种属	Human
特异性	Antibodies were raised in goat to human IgA. There is low background fluorescence in normal tissue in contrast to cases of autoimmune disease where IgA is deposited. No cross reactivity with albumine fractions of human serum. No cross reactivity with albumine fractions of animal species.
经测试应用	适用于: ICC/IF
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	50% glycerol, 0.1 M NaCl, 0.01 M sodium phosphate at pH 7.5.
纯度	Immunogen affinity purified
纯化说明	From the crude polyclonal the cross reactive antibodies were extracted by incubation with Sepharose-bound human IgM and IgG. Specific antibodies were absorbed by incubation with Sepharose-bound human IgA. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized human IgA a minimum of 65% protein bound.
共轭说明	The purified polyclonal was conjugated to fluorescein isothiocyanate isomer I followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
克隆	多克隆
常规说明	F/P ratio is 4-6 mol Fluorescein per 1 mol of goat IgA.

应用

Our [Abpromise guarantee](#) covers the use of **ab8500** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

ICC/IF

应用说明

For immunohistochemical staining of human tissue sections a 1:40 dilution is recommended.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors