

Product datasheet

Goat Anti-Human IgA alpha chain (Biotin) ab97213

概述

产品名称	山羊抗人IgA alpha chain (Biotin)
宿主	Goat
靶标种属	Human
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgA. Cross reactivity with other immunoglobulins and light chains is less than 0.1%.
经测试应用	适用于: IHC-P, ICC/IF, ELISA, WB
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin. Biotinylated antibody was demonstrated with avidin/peroxidase when coated on microtiter wells.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab97213** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/250 - 1/2500.
ICC/IF		1/100 - 1/500.

应用	Ab评论	说明
ELISA		1/10000 - 1/200000. (Primary).
WB		1/10000 - 1/200000. (Chemiluminescent).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors