

## Product datasheet

# Goat Anti-Mouse IgG2a (SureLight® Allophycocyanin) ab72558

### 概述

产品名称	山羊Anti-Mouse IgG2a (SureLight® Allophycocyanin)
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
免疫原	Mouse IgG2a
偶联物	SureLight® Allophycocyanin. Ex: 652nm, Em: 657nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 2mM Sodium Azide Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 100mM Sodium phosphate, pH 7.4
纯度	Size exclusion
纯化说明	Purified using size exclusion chromatography as well as some proprietary methods for purifying the final materials.
共轭说明	Excitation max. $\lambda$ 650 nm. Emission max. $\lambda$ 660 nm.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab72558** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors