

Product datasheet

Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed ab47827

4 References 1 图像

概述

产品名称	山羊抗小鼠IgG+IgM H&L (HRP)预吸附二抗
宿主	Goat
靶标种属	Mouse
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Peroxidase, Mouse IgG and Mouse IgM.
经测试应用	适用于: ICC/IF, ELISA, WB, IHC-P, IHC-Fr, Electron Microscopy
最小交叉反应	Cow, Horse, Human more details
免疫原	Full length native proteins (purified) (Mouse)
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Gentamicin Sulfate Constituents: 10mg/ml BSA (IgG and protease free), 0.15M Sodium chloride, 0.02M Potassium phosphate, pH 7.2
纯度	Immunogen affinity purified
纯化说明	Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption (s) to remove any unwanted reactivities.
克隆	多克隆
同种型	IgG


应用

Our [Abpromise guarantee](#) covers the use of **ab47827** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		1/500 - 1/2500.
ELISA		1/5000 - 1/40000.
WB		1/2000 - 1/10000.
IHC-P		1/500 - 1/2500.
IHC-Fr		1/500 - 1/2500.
Electron Microscopy		Use at an assay dependent dilution.

图片



Ab47827 was used at dilution 1/500 with the primary antibody [ab2869](#) in IHC-P. See the review on [ab2869](#).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat polyclonal Secondary Antibody to Mouse IgG+IgM - H&L (HRP), pre-adsorbed (ab47827)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors