

Product datasheet

Goat Anti-Mouse IgG+IgM H&L (HRP) preadsorbed ab47827

4 References 1 图像

概述

<b>产品名称</b>	山羊抗小鼠IgG+IgM H&L (HRP)预吸附二抗
<b>宿主</b>	Goat
<b>靶标种属</b>	Mouse
<b>特异性</b>	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Peroxidase, Mouse IgG and Mouse IgM.
<b>经测试应用</b>	<b>适用于:</b> ICC/IF, ELISA, WB, IHC-P, IHC-Fr, Electron Microscopy
<b>最小交叉反应</b>	Cow, Horse, Human <a href="#">more details</a>
<b>免疫原</b>	Full length native proteins (purified) (Mouse)
<b>偶联物</b>	HRP

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	pH: 7.20 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
<b>纯度</b>	IgG and protease free
<b>纯化说明</b>	Immunogen affinity purified Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption (s) to remove any unwanted reactivities.
<b>克隆</b>	多克隆
<b>同种型</b>	IgG


应用

Our [Abpromise guarantee](#) covers the use of **ab47827** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		1/500 - 1/2500.
ELISA		1/5000 - 1/40000.
WB		1/2000 - 1/10000.
IHC-P		1/500 - 1/2500.
IHC-Fr		1/500 - 1/2500.
Electron Microscopy		Use at an assay dependent dilution.

## 图片



Ab47827 was used at dilution 1/500 with the primary antibody [ab2869](#) in IHC-P. See the review on [ab2869](#).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat polyclonal Secondary Antibody to Mouse IgG+IgM - H&L (HRP), pre-adsorbed (ab47827)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors