

Product datasheet

Anti-Guanylyl Cyclase beta 1 antibody ab117424

2 图像

概述

产品名称	Anti-Guanylyl Cyclase beta 1 抗体
描述	兔多克隆抗体 to Guanylyl Cyclase beta 1
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow, Dog
免疫原	Synthetic peptide: EDFYEDLDRFEENGTSQESR , corresponding to amino acids 189-207 of Human Guanylyl Cyclase beta 1 (NP_000848.1). Run BLAST with Run BLAST with
阳性对照	Human brain, cerebellum and cortex tissues.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituents: 99% PBS, 1% BSA
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab117424** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		Use a concentration of 10 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标

序列相似性 Belongs to the adenylyl cyclase class-4/guanylyl cyclase family. Contains 1 guanylate cyclase domain.


细胞定位 Cytoplasm.

图片



ab117424, 10 ug/ml, staining Guanylyl Cyclase beta 1 in formalin fixed, paraffin embedded Human brain, cerebellum tissue by Immunohistochemistry followed by biotinylated secondary antibody, alkaline phosphatase-streptavidin and chromogen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Guanylyl Cyclase beta 1 antibody (ab117424)



ab117424, 10 ug/ml, staining Guanylyl Cyclase beta 1 in formalin fixed, paraffin embedded Human brain, cortex tissue by Immunohistochemistry followed by biotinylated secondary antibody, alkaline phosphatase-streptavidin and chromogen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Guanylyl Cyclase beta 1 antibody (ab117424)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors