

Product datasheet

Anti-HA tag antibody (Phycoerythrin) ab72564

2 References

概述

| | |
|-------|---|
| 产品名称 | Anti-HA tag抗体(Phycoerythrin) |
| 描述 | 兔多克隆抗体to HA tag (Phycoerythrin) |
| 偶联物 | Phycoerythrin. Ex: 488nm, Em: 575nm |
| 经测试应用 | 适用于: Flow Cyt |
| 免疫原 | HA tag |
| 常规说明 | Excitation maximum:565>498 nm Emission maximum:578 nm |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | Preservative: 2mM Sodium Azide Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 100mM Sodium phosphate, pH 7.4 |
| 纯度 | Size exclusion |
| 纯化说明 | Size exclusion chromatography as well as some proprietary methods for purifying the final materials were used. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab72564** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----------|------|--|
| Flow Cyt | | Use at an assay dependent concentration. ab37409 -Rabbit polyclonal IgG, is suitable for use as an isotype control with this antibody. |

靶标

相关性

Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein. This tag facilitates the detection, isolation, and purification of the proteins.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors