

# Anti-HIV1 gp120 antibody [HIV1gp120] ab13411

### 概述

产品名称	Anti-HIV1 gp120抗体[HIV1gp120]
描述	小鼠单克隆抗体[HIV1gp120] to HIV1 gp120
宿主	Mouse
经测试应用	适用于: ELISA
免疫原	r.gp120 (MN strain).
常规说明	Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. If you would like help finding an alternative product, please do not hesitate to contact our scientific support team.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: PBS, pH 7.4
纯度	Ion Exchange Chromatography
纯化说明	Antibody was purified by ion exchange chromatography.
克隆	单克隆
克隆编号	HIV1gp120
同种型	IgG1

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用 ab13411 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

## 靶标

---

相关性	HIV1 is equipped with the envelope gp160 glycoprotein for interaction with Langerhans cells (LCs) and dendritic cells (DCs), the members of the innate immune system, which confront the virus at the portal of virus entry in the human body. These cells are equipped with receptors by which they bind and endocytose the virus. The gp120 glycoprotein is used for binding to CD4 receptor and CCR5 co-receptor of T helper 2 (Th2) cells, and is able to induce FcepsilonRI(+) hematopoietic cells to produce IL4, which inactivates the host adaptive immune response.
细胞定位	Envelope surface glycoprotein.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors