

Product datasheet

Anti-HIV1 tat antibody (FITC) ab43016

1 图像

概述

产品名称	Anti-HIV1 tat抗体(FITC)
描述	兔多克隆抗体to HIV1 tat (FITC)
宿主	Rabbit
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	This antibody is specific for Tat protein and GST-Tat fusion protein.
经测试应用	<b>适用于:</b> Flow Cyt
免疫原	Recombinant full length protein expressed in E. coli.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: None Constituents: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab43016** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Flow Cyt

应用说明	Flow Cyt: 1/200. For general use, dilute at 1:1 (V/V) in PBS.  Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
------	---

## 靶标


### 相关性

HIV1 tat (Transactivator of Transcription) protein is a pleiotropic factor that induces a broad range of biological effects in numerous cell types. At the HIV promoter, tat is a powerful transactivator of gene expression which acts by inducing chromatin remodeling and by recruiting elongation-competent transcriptional complexes on to the viral LTR.

### 细胞定位

Nuclear. Upon stimulation with EDN1, it is exported from the nucleus to the perinuclear region.

## 图片



ab43016 at a 1/200 dilution.

Flow Cytometry - Anti-HIV1 tat antibody (FITC)  
(ab43016)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors