


Product datasheet

Anti-HSPC302 antibody ab96188

1 图像

概述

产品名称	Anti-HSPC302抗体
描述	兔多克隆抗体to HSPC302
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Zebrafish 
免疫原	Recombinant fragment containing a sequence corresponding to a region within amino acids 394-601 of Human HSPC302 (NP_149106).
阳性对照	HepG2 and Molt4 whole cell lysates; A431, H1299, HeLa and Raji lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab96188** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 101 kDa.

靶标

Anti-HSPC302 antibody 图像



All lanes : Anti-HSPC302 antibody
(ab96188) at 1/1000 dilution

Lane 1 : HepG2 whole cell lysate

Lane 2 : Molt4 whole cell lysate

Lysates/proteins at 30 µg per lane.

Predicted band size : 101 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors