

## Product datasheet

# Anti-Haemophilus influenzae a antibody ab79394

### 概述

产品名称	Anti-Haemophilus influenzae a抗体
描述	兔多克隆抗体to Haemophilus influenzae a
宿主	Rabbit
特异性	No cross reactivity with other serotypes observed.
经测试应用	适用于: Agglutination
种属反应性	Haemophilus influenzae a
免疫原	Tissue/ cell preparation of Haemophilus influenzae a (heat inactivated)
阳性对照	Antigenic suspension.
常规说明	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.08% Sodium Azide Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab79394** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Agglutination		Use at an assay dependent concentration. See attached protocol.

## 靶标

---

**相关性** Haemophilus influenzae is a small, non-motile Gram negative bacterium found naturally in the nasopharynx of approximately 75 percent of healthy children and adults. First described by Pfeiffer in 1892, some forty years later Pittman identified six capsular polysaccharides (types) of H. influenzae (a, b, c, d, e, and f). Similar to Haemophilus influenzae b, most cases of systemic Haemophilus influenzae a infection are in young children with meningitis being the most common disease.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors