

Product datasheet

# Anti-Hantavirus Puumala strain antibody [1D11] ab116453

概述

产品名称	Anti-Hantavirus Puumala strain抗体[1D11]
描述	小鼠单克隆抗体[1D11] to Hantavirus Puumala strain
特异性	ab116453 is specific to Hantavirus Puumala strain.
经测试应用	适用于: WB, ELISA
免疫原	Recombinant Hantavirus Puumala strain.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 99% PBS
纯度	Protein A purified
纯化说明	Purity is > 95% by SDS-PAGE.
克隆	单克隆
克隆编号	1D11
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab116453** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 1 kDa.
ELISA		Use at an assay dependent dilution.

靶标

## 相关性

Hantaviruses are rodent-borne agents that cause severe human diseases. The viruses constitute a group of antigenetically and genetically related envelope viruses with a tripartite, negative-sense RNA genome. In northern, western and central Europe the Puumala virus (PUUV) is a predominant serotype. It is the etiological agent of nephropathia epidemica, a mild form of hemorrhagic fever with renal syndrome.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors