

## Product datasheet

# Human CITED2 peptide ab22995

### 概述

产品名称 人CITED2多肽

### 描述

性质 Synthetic

### 氨基酸序列

种属 Human

序列 DFVCKQQPSRVSC

氨基酸 258 to 270

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab22995** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking

形式 Liquid

### 制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### 常规信息

功能 Interferes with the binding of transcription factors HIF-1a and STAT2 to p300/CBP.

序列相似性 Belongs to the CITED family.

细胞定位 Nucleus.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors