

Product datasheet

Human FAM38B peptide ab93376

概述

产品名称 人FAM38B多肽

描述

性质 Synthetic

氨基酸序列

种属 Human

序列 MIAGNSTESS KTPVTIEKIY PYYVKAPSDS
NSKPIKQLLS ENNFMDITII

氨基酸 360 to 409

技术指标

Our [Abpromise guarantee](#) covers the use of **ab93376** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

Preservative: None

Constituents: ddH2O

常规信息

功能 Component of a mechanosensitive channel required for rapidly adapting mechanically activated (MA) currents. Required for Merkel-cell mechanotransduction. Plays a major role in light-touch mechanosensation.

疾病相关 Arthrogryposis, distal, 5
Arthrogryposis, distal, 3
Marden-Walker syndrome

序列相似性	Belongs to the PIEZO (TC 1.A.75) family.
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors