

# Human FOXN1 peptide ab45510

## 概述

产品名称 人FOXN1肽

## 描述

性质 Synthetic  
 氨基酸序列  
 种属 Human  
 序列 C-SVYLSPSSKPVALA  
 氨基酸 635 to 648

## 技术指标

Our [Abpromise guarantee](#) covers the use of **ab45510** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking - Blocking peptide for Anti-FOXN1 antibody ([ab26069](#))  
 形式 Liquid

## 制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## 常规信息

功能 Transcriptional regulator involved in development.  
 组织特异性 Expressed in thymus.  
 疾病相关 Defects in FOXN1 are the cause of T-cell immunodeficiency congenital alopecia and nail dystrophy (TIDAND) [MIM:601705]. A disorder characterized by the association of congenital alopecia, severe T-cell immunodeficiency, and ridging and pitting of all nails.  
 序列相似性 Contains 1 fork-head DNA-binding domain.  
 细胞定位 Nucleus.

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors