

Product datasheet

Human KF1 peptide ab23039

概述

产品名称 人KF1多肽

描述

性质 Synthetic

氨基酸序列

种属 Human

序列 C-YAQHQPLSNDVPS

氨基酸 673 to 685

技术指标

Our [Abpromise guarantee](#) covers the use of **ab23039** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

常规信息

功能 Acts as an E2-dependent E3 ubiquitin-protein ligase, probably involved in the ER-associated protein degradation pathway.

组织特异性 Highly expressed in the normal cerebellum but not in the cerebral cortex.

通路 Protein modification; protein ubiquitination.

序列相似性 Contains 1 RING-type zinc finger.

细胞定位 Endoplasmic reticulum membrane.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors