

## Product datasheet

# Human Substance P peptide ab42295

### 概述

**产品名称** 人Substance P肽

### 描述

**性质** Synthetic

### 氨基酸序列

**种属** Human

**序列** RPRPQQFFGLM

**氨基酸** 58 to 68

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab42295** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**应用** Blocking - Blocking peptide for Anti-Substance P antibody [SP-DE4-21] ([ab14184](#))

**形式** Lyophilised

### 制备和贮存

**稳定性和存储** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Preservative: None

Constituents: 0.1% Acetic acid

**复溶** Resuspend in 0.1% acetic acid for required concentration

### 常规信息

**功能** Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodilators and secretagogues, and contract (directly or indirectly) many smooth muscles.

**序列相似性** Belongs to the tachykinin family.

**细胞定位** Secreted.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors