

Product datasheet

Human Syntaxin 6 peptide ab48062

概述

产品名称 人Syntaxin 6肽

描述

性质 Synthetic

氨基酸序列

种属 Human

序列 C-QALAERKNRQA

氨基酸 112 to 122

技术指标

Our [Abpromise guarantee](#) covers the use of **ab48062** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking - Blocking peptide for Anti-Syntaxin 6 antibody ([ab48016](#))

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

常规信息

功能 Involved in intracellular vesicle trafficking.

序列相似性 Belongs to the syntaxin family.  
Contains 1 t-SNARE coiled-coil homology domain.

翻译后修饰 Phosphorylated upon DNA damage, probably by ATM or ATR.

细胞定位 Golgi apparatus membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors