

## Product datasheet

# Human TXNDC peptide ab45442

### 概述

**产品名称** 人TXNDC肽

### 描述

**性质** Synthetic

### 氨基酸序列

**种属** Human

**序列** RSLGPSLATDKS

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab45442** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**应用** Blocking - Blocking peptide for Anti-TXNDC antibody ([ab37876](#))

**形式** Liquid

### 制备和贮存

**稳定性和存储** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Preservative: None

### 常规信息

**功能** May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions.

**组织特异性** Ubiquitous. Highly expressed in kidney, liver, placenta and lung.

**序列相似性** Contains 1 thioredoxin domain.

**细胞定位** Membrane. Endoplasmic reticulum membrane. Predominantly found in the endoplasmic reticulum.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors