

### Anti-ITPR2 antibody ab5805

[4 References](#) [2 图像](#)

#### 概述

<b>产品名称</b>	Anti-ITPR2抗体
<b>描述</b>	兔多克隆抗体to ITPR2
<b>宿主</b>	Rabbit
<b>经测试应用</b>	<b>适用于:</b> ICC, IHC-Fr <b>不适用于:</b> WB
<b>种属反应性</b>	<b>与反应:</b> Rat, Human <b>预测可用于:</b> Cow 
<b>免疫原</b>	Synthetic peptide corresponding to Rat ITPR2 aa 317-334. Sequence: ELNPDYRDAQNEGKTVRD  (Peptide available as <a href="#">ab5842</a> )
<b>阳性对照</b>	IHC-Fr: Rat brain tissue. ICC: SH-SY5Y cells.
<b>常规说明</b>	The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.  If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As


 [Run BLAST with](#)

 [Run BLAST with](#)

#### 性能


<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>存储溶液</b>	Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 99% PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆





Immunocytochemistry - Anti-ITPR2 antibody (ab5805)

4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilised SH-SY5Y (human neuroblastoma cell line from bone marrow) cells labeling ITPR2 using ab5805 at 2 µg/ml in 0.1% BSA for 3 hours at RT, followed by a Goat Anti-Rabbit (H&L) AlexaFluor® secondary antibody (Green (a)). The nuclear counter stain is DAPI (Blue (b)). F-actin is stained with AlexaFluor® 555-Rhodamine (Red (c)). Panel D shows the merged image. Panel (e) is a secondary antibody only control.



Immunohistochemistry (Frozen sections) - Anti-ITPR2 antibody (ab5805)

Shows immunolocalization of ITPR2 in rat brain using ab5805.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors