

## Product datasheet

# Anti-Interferon alpha 2 antibody [N27] ab24099

## 2 References

### 概述

产品名称	Anti-Interferon alpha 2抗体[N27]
描述	小鼠单克隆抗体[N27] to Interferon alpha 2
宿主	Mouse
特异性	Specific for Interferon alpha 2 and does not crossreact with Interferon alpha 1.
经测试应用	<b>适用于:</b> ELISA, Blocking, WB <b>不适用于:</b> IHC-P
种属反应性	<b>与反应:</b> Human
免疫原	E.coli derived recombinant full length Interferon alpha 2c (Human)
表位	ab24099 reacts with linear peptide 43aa to 53aa of the 172aa total for Interferon alpha 2.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium azide Constituent: PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	N27
骨髓瘤	NS1
同种型	IgG1
轻链类型	kappa

### 应用

Our [Abpromise guarantee](#) covers the use of **ab24099** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
Blocking		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 22 kDa.

**应用说明** Is unsuitable for IHC-P.

## 靶标

**功能** Produced by macrophages, IFN-alpha have antiviral activities.

**序列相似性** Belongs to the alpha/beta interferon family.

**细胞定位** Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors