

## Product datasheet

# Anti-KCNK2 antibody ab90855

★★★★☆ 4 Abreviews 2 References 1 图像

### 概述

产品名称	Anti-KCNK2抗体
描述	兔多克隆抗体to KCNK2
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow
免疫原	Synthetic peptide within the N terminal region (amino acids 1-100) of Human KCNK2 (UniProt ID: Q9NRT2). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Human brain membrane lysate

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab90855** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

**应用说明** WB: Use at a concentration of 2 µg/ml. Predicted molecular weight: 47 kDa. It is essential to run against membrane preps.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标

功能	Outward rectifying potassium channel.
序列相似性	Belongs to the two pore domain potassium channel (TC 1.A.1.8) family.
细胞定位	Membrane.

## 图片



Anti-KCNK2 antibody (ab90855) at 2  $\mu$ g/ml +  
Human brain membrane lysate

**Predicted band size : 47 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors