

Product datasheet

Anti-KLHDC7B antibody ab121102

2 图像

概述

产品名称	Anti-KLHDC7B抗体
描述	兔多克隆抗体to KLHDC7B
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to amino acids 348-489 of Human KLHDC7B.
阳性对照	RT4, U251 MG, A431 cell lysates; Human urinary bladder tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121102** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标

### 序列相似性

Contains 5 Kelch repeats.

## 图片



**All lanes :** Anti-KLHDC7B antibody (ab121102) at 1/250 dilution

**Lane 1 :** RT4 cell lysate

**Lane 2 :** U251 MG cell lysate


**Lane 3 :** A431 cell lysate

**Lane 4 :** Human Liver cell lysate

**Lane 5 :** Human Tonsil cell lysate

Developed using the ECL technique.

10-20% SDS PAGE



ab121102 at 1/100 dilution staining KLHDC7B in Paraffin-embedded Human urinary bladder tissue by Immunohistochemistry.

Note: moderate membranous positivity in urothelial cells.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors