

Product datasheet

Anti-LCMV antibody [M104] ab31774

1 Abreviews

概述

产品名称	Anti-LCMV抗体[M104]
描述	小鼠单克隆抗体[M104] to LCMV
宿主	Mouse
特异性	This antibody reacts with LCMV infected cells (e.g. Vero, HeLa). The originator of clone M104 reported that this clone recognizes the nucleocapsid protein of LCMV.
经测试应用	适用于: ICC/IF, WB, IHC-Fr
免疫原	Full length protein: LCMV isolated from human HeLa cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Preservative: 0.09% Sodium Azide. Constituents: Tissue culture supernatant.
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	M104
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab31774** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

靶标

相关性

Lymphocytic choriomeningitis virus (LCMV), the first member of the arenavirus family to be isolated, is the causative agent of a zoonosis acquired from chronically viremic mice or hamsters. LCMV primarily infects wild mice and it is estimated that 5% of all wild mice carry LCMV. Mice can asymptotically carry and shed the virus in saliva, urine and faeces. People and other animals become infected through contact with these secretions or by inhalation of dried particles from them. The clinical spectrum of acquired human LCMV infection ranges from inapparent and asymptomatic to, in rare instances, severely symptomatic, systemic, and fatal central nervous system (CNS) disease. Intrauterine LCMV infection has resulted in fetal or neonatal death, as well as hydrocephalus and chorioretinitis in infants.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors