

Product datasheet

Anti-LOC285331 antibody ab118801

1 图像

概述

产品名称	Anti-LOC285331抗体
描述	兔多克隆抗体to LOC285331
宿主	Rabbit
经测试应用	适用于: WB 不适用于: IP
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to a region within amino acids 898-948 of Human LOC285331 (NP_001135419.1).
阳性对照	293T and Jurkat whole cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate Note: pH 7 to 8
纯度	Immunogen affinity purified
纯化说明	ab118801 was affinity purified using an epitope specific to LOC285331 immobilized on solid support.
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab118801** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000 - 1/10000. Predicted molecular weight: 109 kDa.

靶标

图片



Western blot - Anti-LOC285331 antibody (ab118801)

All lanes : Anti-LOC285331 antibody
(ab118801) at 0.1 $\mu\text{g/ml}$

Lane 1 : 293T whole cell lysate at 50 μg

Lane 2 : 293T whole cell lysate at 15 μg

Lane 3 : Jurkat whole cell lysate at 50 μg

Developed using the ECL technique.

Predicted band size: 109 kDa

Exposure time: 3 minutes

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors