

Product datasheet

Anti-LOC653186 antibody ab69849

1 图像

概述

产品名称	Anti-LOC653186抗体
描述	兔多克隆抗体to LOC653186
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Synthetic peptide: MKVINAKLDG FYSFPIFLFQ FLQATAQEEG IFECADPKLA ISAWTFRDL, corresponding to amino acids 1-50 of Human LOC653186 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	HepG2 cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab69849** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use a concentration of 2.5 µg/ml. Detects a band of approximately 10 kDa (predicted molecular weight: 10 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


应用	Ab 评论	说明
----	-------	----

ELISA Use at an assay dependent concentration.  
ELISA titre using peptide based assay: 1/312500.

## 靶标

**相关性** The function of LOC653186 remains unknown.

## 图片



Anti-LOC653186 antibody (ab69849) at 2.5 µg/ml + HepG2 cell lysate at 10 µg/ml

### Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

**Predicted band size** : 10 kDa

**Observed band size** : 10 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors