

Product datasheet

Anti-LPPR3 antibody ab121021


概述

产品名称

Anti-LPPR3抗体

描述

山羊多克隆抗体to LPPR3

 该产品是 **快线抗体**，通过亲和纯化得到。对抗原肽有高灵敏度及特异性 (ELISA测试)。

[Read the terms of use »](#)

宿主

Goat

种属反应性

预测可用于: Rat, Human 

免疫原

Synthetic peptide:

C-TQRGHDSVYQQN

, corresponding to internal sequence amino acids 335-346 of Human LPPR3 (NP_079164.1).

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml).

性能

形式

Liquid

存放说明

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

存储溶液

pH: 7.30

Preservative: 0.02% Sodium azide

Constituents: 99% Tris buffered saline, 0.5% BSA

纯度

Immunogen affinity purified

纯化说明

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

克隆

多克隆

同种型

IgG

应用

Fast track antibodies constitute a diverse group of products that have been released to accelerate your research, but are not yet fully characterized. They have all been affinity purified and show high titre values against the immunizing peptide (by ELISA).

[Fast track terms of use](#)

应用	Ab 评论	说明
----	-------	----

WB		Use at an assay dependent concentration. Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
----	--	--

靶标

序列相似性

Belongs to the PA-phosphatase related phosphoesterase family.

细胞定位

Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors