

Product datasheet

Anti-LRRC19 antibody ab106657

2 图像

概述

产品名称	Anti-LRRC19抗体
描述	兔多克隆抗体to LRRC19
宿主	Rabbit
特异性	ab106657 is predicted not to cross-react with other leucine rich repeat family proteins or TLRs.
经测试应用	适用于: ICC, ICC/IF, WB
种属反应性	与反应: Mouse, Human
免疫原	A 17 amino acid synthetic peptide near the C terminus of Human LRRC19 (NP_075052).
阳性对照	K562 cells and cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106657** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC		Use a concentration of 2.5 µg/ml.
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 42 kDa.

靶标


序列相似性

Contains 5 LRR (leucine-rich) repeats.
Contains 1 LRRCT domain.

细胞定位


Membrane.

图片



Anti-LRRC19 antibody (ab106657) at 1 μ g/ml
+ K562 cell lysate at 15 μ g

Predicted band size: 42 kDa



ab106657, at 2.5 μ g/ml, staining LRRC19 in K562 cells by Immunocytochemistry.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors