

Product datasheet

Anti-LRRC56 antibody ab122506

2 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-LRRC56抗体 |
| 描述 | 兔多克隆抗体to LRRC56 |
| 经测试应用 | 适用于: IHC-P, ICC/IF |
| 种属反应性 | 与反应: Human |
| 免疫原 | antigen, corresponding to amino acids 29-133 of Human LRRC56 (Q8YIG6). |
| 阳性对照 | Human testis tissue. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at -20°C. |
| 存储溶液 | pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab122506** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| 应用 | Ab评论 | 说明 |
|--------|------|---|
| IHC-P | | 1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |
| ICC/IF | | Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment |

靶标

序列相似性


Contains 5 LRR (leucine-rich) repeats.
Contains 1 LRRCT domain.

图片



Immunofluorescent staining of Human cell line U-251 MG shows positivity in nucleus but not nucleoli. Recommended concentration of ab122506 1-4 $\mu\text{g/ml}$. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-LRRC56 antibody (ab122506)



ab122506, at 1/200 dilution, staining LRRC56 in formalin-fixed, paraffin-embedded Human testis tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-LRRC56 antibody (ab122506)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors