

Product datasheet

Anti-LUC7L2 antibody ab104004

2 图像

概述

产品名称	Anti-LUC7L2抗体
描述	兔多克隆抗体to LUC7L2
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 348-377 of Human LUC7L2.
阳性对照	MCF-7 cell lysate; Human brain tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	ab104004 was purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab104004** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 47 kDa.

应用	Ab评论	说明
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标


功能	May bind to RNA via its Arg/Ser-rich domain.
序列相似性	Belongs to the Luc7 family.
细胞定位	Nucleus speckle. Nucleus, nucleoplasm. Colocalizes with SCNM1 and SNRNP70 in nuclear speckles.
形式	Cellular localisation: nucleus

## 图片



Anti-LUC7L2 antibody (ab104004) at 1/100 dilution + MCF-7 cell lysate at 35 µg

**Predicted band size:** 47 kDa



ab104004, at 1/50 dilution, staining LUC7L2 in formalin-fixed, paraffin-embedded Human brain tissue by Immunohistochemistry using a peroxidase conjugated secondary antibody and DAB staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors