

Product datasheet

M-CSF overexpression 293T lysate (whole cell) ab94076

2 图像

概述

产品名称	M-CSF overexpression 293T裂解液(全细胞)
描述	M-CSF overexpression裂解液
常规说明	ab94076 is a 293T cell transfected lysate in which Human MCSF has been transiently over-expressed using a pCMV-MCSF plasmid. The lysate is provided in 1X Sample Buffer. <b>Note:</b> For more details about how the transfected lysate was prepared <a href="#">view preparation notes</a>
经测试应用	适用于: WB

性能

形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 6.80 Constituent: 100% 1x Sample Buffer
纯度	Whole Cell Lysate
背景	Function: Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. CSF-1 induces cells of the monocyte/macrophage lineage. It plays a role in immunological defenses, bone metabolism, lipoproteins clearance, fertility and pregnancy. PTM: Glycosylation and proteolytic cleavage yield different soluble forms. A high molecular weight soluble form is a proteoglycan containing chondroitin sulfate. Isoform 1 is N- and O-glycosylated. Isoform 3 is N-glycosylated.

应用

Our [Abpromise guarantee](#) covers the use of **ab94076** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



SDS-PAGE - MCSF 293T Transfected Lysate -  
(positive control) (ab94076)

ab94076 at 15µg/lane on an SDS-PAGE gel.



Western blot - MCSF 293T Transfected Lysate -  
(positive control) (ab94076)

**All lanes** : Anti-MCSF antibody (ab66236) at  
1/500 dilution

**Lane 1** : M-CSF overexpression 293T lysate  
(whole cell) (ab94076)

**Lane 2** : 293T non transfected lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

Goat Anti-mouse IgG (H and L) HRP  
conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors