

Product datasheet

Anti-MARCH4 antibody ab78005

1 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-MARCH4抗体 |
| 描述 | 兔多克隆抗体to MARCH4 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Mouse, Rat, Horse, Chicken, Dog, Human, Chimpanzee, Zebrafish |
| 免疫原 | Synthetic peptide, corresponding to a portion of the amino acids 300-350 of Human MARCH4 |
| 阳性对照 | Liver tissue. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle. |
| 存储溶液 | Preservative: 0.05% Sodium Azide Constituents: 0.2% Gelatin, PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab78005** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|---|
| WB | | Use a concentration of 0.5 - 2 µg/ml. Predicted molecular weight: 46 kDa. |
| IHC-P | | Use a concentration of 5 µg/ml. |

靶标

| | |
|-------------|---|
| 相关性 | MARCH4, an E3 ubiquitin protein ligase, may mediate ubiquitination of MHC I and CD4, and promote their subsequent endocytosis and sorting to lysosomes via multivesicular bodies. E3 ubiquitin protein ligases accept ubiquitin from an E2 ubiquitin conjugating enzyme in the form of a thioester and then directly transfer the ubiquitin to targeted substrates. |
| 细胞定位 | Golgi apparatus membrane; Multipass membrane protein. |

图片



ab78005, at 5µg/ml, staining MARCH4 in liver tissue by Immunohistochemistry using formalin-fixed, paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MARCH4 antibody (ab78005)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors